CoaCPN (English) <div>Control Plasma Normal</div> Reagents suitable for CoaData and CoaLAB analyzer	 
--	---

	210-24-050-00		8x1mL
---	----------------------	---	--------------

INTENDED PURPOSE

CONTROL PLASMA NORMAL is a material that guarantees the adequacy of the performances and the procedure of use of the kit for the quantitative determination with specific methods of coagulation parameters in plasma for the reagents:

- CoaPT Prothrombin time  **210-24-010-00**
- CoaPTT Activated Partial Thromboplastin Time  **210-24-020-00**
- CoaFIB-Fibrinogen Kit (Claus Method)  **210-24-040-00**


PRINCIPLE

To use it follow the procedures shown on the Information for Use of the kits for the determination of each single parameter

PRECAUTION FOR USE

1.This product has been formulated for in vitro diagnostic use.
2.A proportional variation of the reaction volumes does not change the result.
3.DO NOT mix Reagents from different Production lots.
4.In addition to the possible risk indications, the Reagent can contain preservatives, which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However, it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes, and skin; to use reagents according to good laboratory practice.
On the Safety Data Sheet are detailed the operating procedures for the manipulation of this product. Safety Data Sheet should be supplied on request.

ATTENTION!

A) The reagent can be used with manual, mechanical, photometric and nephelometric clot detection systems. The automated determinations must be performed according to specific instructions attached to the instrument used.
B) This product comes from donors and has to be used as a potential transfer of infective pathologies. This product has been tested and found to be negative for HIV, HCV and HBsAg antibodies by an approved method. However, as there is no tested method that can offer complete assurance that HIV I-II, Hepatitis virus and other infectious agents are absent, the **210-24-050-00** and all samples have to be handled as potential transmitter of infectious diseases, using the procedures detailed on Safety Data Sheets.
C) The reagent must be used **ONLY** for the intended purpose, by expert and trained people and according to good laboratory practice.
D) The clinical diagnosis cannot be done correctly using the result of only one test but have to be done integrating critically the results of different laboratory tests and clinical data.
E) A series of factors, such as ambient temperature, the temperature of the working reagents, the accuracy of the washings, the type of coagulometer and the distilled water characteristics, can affect the performance of the test.
F) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.
G) All the precautions normally used in the laboratory must be respected for reagents handling.
All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

REAGENTS kit composition:	 210-24-050-00
R1 - 	8 x 1 mL
Human stabilized plasma	
(See the values table for each parameter)	

STABILITY

The Reagents are stable up to the expiry date mentioned on the labels, if closed and stored at 2-8°C in their intact primary container, if not exposed to heat sources and/or pressure variations.
In case of damaging of the primary container organize the waste disposal.

PREPARATION OF THE WORKING REAGENT

Open one vial very carefully, avoiding to lose the lyo material. Add exactly 1 mL of distilled water; close bottle carefully and let stand for 30 minutes. Dissolve contents completely by swirling gently, avoiding the foam formation. Mix kindly before use and let the reagent reach the room temperature before use. Controls should be used following the procedures required for the patient sample. Refer to the instructions for use supplied in the kit of each single parameter and to the instruction manual of the coagulometer. Close immediately after handling. The Reagents have to be used correctly, to avoid contamination. An incompetent handling relieves us from any responsibility.

STABILITY AFTER RECONSTITUTION

The product is stable:
- 24 hours at 2-8°C for PT, aPTT, Fibrinogen and Antithrombin.
- 4 hours open vial in use +18 - 26°C
- 1 month at -20°C for PT, aPTT, Fibrinogen and Antithrombin.
After reconstitution FREEZE only ONE TIME, DO NOT REPEAT FREEZING.
In case of microbial contamination, the preparations must be eliminated.

MATERIAL REQUIRED BUT NOT PROVIDED

Standard laboratory equipment.
Micropipettes to deliver from 3 to 1000 µL.
Disposable micropipettes tips.
Plastic test tubes for sample dilution.
Coagulation tubes.
Stopwatch or timer.
Water bath 37°C.
Distilled water, Reagents and Calibrators.
Coagulometers.

ANALYTICAL PROCEDURE ON MANUAL and AUTOMATIC COAGULOMETER

Let reagents reach the working temperature before use. Use controls as natural wild human plasma sample.

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ATTENTION!

The kit is tested on Coagulometer CoaLAB 1000, CoaData 4004 (both LABiTec) and ACL 7000 (I.L.). Applications on automatic / semi-automatic or manual coagulometer from other manufacturers may be totally different from what we experienced.

CALCULATION OF RESULTS

Please refer to the IFU of individual coagulation reagents.

ANALYTICAL PERFORMANCES

(validate on CoaLAB 1000, CoaData 4004 and ACL 7000 (I.L.)
The analytical performances of parameters are detailed on the relative Information for Use.

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. **LABiTec** disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of **LABiTec**'s products for their specific applications. In no event will **LABiTec** be liable for any consequential damages arising out of aforesaid express warranty.

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


REFERENCES

See References in IFU of each reagent.

<p>Symbols key / Symbolschlüssel / Interpretazione simboli / Clave de los Símboles</p>	
	CE Mark / CE Symbol / Marchio CE / Marca CE
	In Vitro Diagnostic Medical Device / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Catalogue Number / Bestellnummer / Numero di catalogo / Número de catálogo
	Refer to instruction for use / Bedienungsanleitung beachten / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Biological Risk / Biologisches Risiko / Rischio biologico / Riesgo biológico
	Manufactured by / Hersteller / Fabricante / Fabricante
	Distributor/Händler/Distributore/Distribuidora
	Temperature limitations / Temperaturbegrenzung / Limiti di temperature / Limite de temperatura
	Contents sufficient to "n" tests / Inhalt ausreichend für „n“Tests / Contenuto sufficiente per numero test / Contenidos suficientes para "n" pruebas
	Lot Number / Chargennummer / Codice del lotto / Código de lote
	Expiration Date (YYYY/MM) / Verwendbar bis / Utilizare entro / Fecha de caducidad
	Contents / Inhalt / Contenuto / Contenido
	Control / Kontrolle / Controllo / Control
	Calibrator / Kalibrator / Calibratore / Calibrador
	Reconstitution Volume / Rekonstitutionsvolumen / Volume di ricostituzione / Volumen de reconstitución
	Respiratory sensitisation, hazard categories 1, 1A, 1B – Germ cellmutagenicity, hazrard 1A, 1B, 2 – Carcinogenicity, hazard categories 1A, 1B, 2 – Reproductive toxicity, hazard 1A, 1B, 2 – Specific Target OrganToxicity – Singel exposure, hazard categories 1, 2 – SpecificTarget Organ Toxicity – Repeated exposure, hazard categories 1, 2 Aspiration hazard, hazard category 1 (Health Hazard)
	Acute toxicity (oral, dermal, inhalation), hazard category 4 - Skin irritation, hazard category 2. Eye irritation, hazard category 2 - Skin sensitisation, hazard categories 1, 1A - Specific Target Organ Toxicity — Single exposure, hazard category 3- Respiratory tract irritation. Narcotic effects (Health Hazard)

	BEN – BIOCHEMICAL ENTERPRISE S.r.l. - via Toselli, 4 - 20127 Milano – Italy Tel: +39 022871829/93 – Fax: +39 022890853 – www.bensrl.it – info@bensrl.it
---	---

Rev. 210-24-050-00_04.24

	BEN – BIOCHEMICAL ENTERPRISE S.r.l. - via Toselli, 4 - 20127 Milano – Italy Tel: +39 022871829/93 – Fax: +39 022890853 – www.bensrl.it – info@bensrl.it
---	---

	
	LABiTec GmbH An der Strusbek 6, 22926 Ahrensburg, Germany www.labitec.com, T: +49410247950, info@labitec.de