CoaPTT (English)

Activated Partial Thromboplastin Time
Reagents suitable for CoaDATA and CoaLAB analyzer

REF 210-24-020-00 **CONT** 8x4mL

INTENDED PURPOSE (see References 1, 4, 7)

The Activated Partial Thromboplastin Time (aPTT) is the most important method for monitoring the "intrinsic" pathway of coagulation of the blood and anticoagulant therapy with heparin. The ability of blood to form a fibrin clot by means of the intrinsic pathway requires phospholipids, calcium and a contact activator negatively charged, ellagic acid. So aPTT is sensitive to the decrease in the concentration of coagulation factors of the intrinsic and common, of the contact phase and of the anticoagulant effects of heparin. A prolongation of aPTT is generally related to the decrease of one or more Factors (deficiency of XII, XI, X, IX, V, II and fibrinogen); by functional deficits (deficiency of vitamin K, liver disease); the effect of an anticoagulant such as heparin; by the presence of an inhibitor.

PRINCIPLE

The CoaPTT reagent consists of stabilized and optimized rabbit brain phospholipids as a platelet substitute and a soluble plasma activator, ellagic acid, for optimal activation of the contact phase of coagulation. The sample plasma is placed in contact with the reagent containing an optimal number of phospholipids and negatively charged contact activator; incubation at 37°C for a defined period allows activation of the intrinsic coagulation pathway. The addition of calcium ions to the reaction mixture initiates coagulation and the time required for fibrin clot formation is determined. The optimization of the CoaPTT Reagent allows to obtain better or comparable performances to those of the market leading Reagents that use ellagic acid: regardless of whether they use selected or synthetic phospholipids or soy phospholipids.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. In addition to the possible risk indications, the Reagent can contain preservatives, which total concentration is lower than the limits mentioned in Dir 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However, it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

- A) The Reagent can be used with manual, mechanical, photometric and nephelometric clot detection systems. The automated determinations must be performed according to specific instructions attached to the instrument used.
 B) Very deep attention must be given to interfering substances: certain drugs and other substances may influence levels of aPTT or aPTT assay (see References 2).
 C) The Reagent must be used ONLY for the intended purpose, by expert and trained people and according to good laboratory practice.
 D) A clinical diagnosis cannot be correctly made based on the result of a single test, but must be made by critically integrating the results of different laboratory tests with different clinical data.
 E) A series of factors, such as ambient temperature, the temperature of the working reagents, the accuracy of the washings, the type of coagulometer and the distilled water characteristics, can affect the performances of the test.
 F) For the handling of reagents, observe the precautions normally taken in the laboratory. All the calibrators and controls are human samples, so potentially infectious; must therefore be taken all appropriate protective measures in order to avoid any potential biohazard.

REAGENTS
Components of the kit: **REF** 210-24-020-00
APTT Reagent, liquid 8 x 4 mL

Rabbit brain cephalin
Kaolin
Ellagic acid
Excipients and stabilizers

STORAGE AND STABILITY

The CoaPTT Reagent is stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container, if not exposed to heat sources and/or pressure variations. In case the primary container is damaged, it should be disposed of properly. The CoaPTT is stable according to the following specification.

Unopened	2 - 8°C	Refer to vial label
Open Vial (in original vial)	2 - 8°C	30 Days
On-board CoaLAB 1000	18 - 26°C	5 days

At the end of the working cycles is suggested to store the reagent to 2-8°C in the original vial for better stability.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right test performances use following kits (see the relative instructions for use (IFU)):

CoaCAL COAGULATION CALIBRATOR **REF** 210-24-070-00
CoaCPN CONTROL PLASMA NORMAL **REF** 210-24-050-00
CoaCPA CONTROL PLASMA ABNORMAL **REF** 210-24-060-00

Each laboratory should establish its own range of control that determines the acceptable variation of the day-to-day performance for each control plasma. The calibration curve must be repeated for each new lot of reagents and / or calibrator.

PREPARATION OF THE WORKING REAGENT

The CoaPTT Reagent is supplied ready to use. Before use, the container must be mixed several times by inversion to ensure the homogeneity of the reagent. A lipid/ellagic acid sediment can form if the containers are left standing for a long time. Mix gently before use: DO NOT shake. NEVER freeze the CoaPTT Reagent. Close immediately after use. Reagents must be used correctly to avoid contamination. Improper handling releases us from any responsibility.

MATERIALS REQUIRED BUT NOT PROVIDED

CoaCC CALCIUM CHLORIDE, liquid **REF** 210-10-030-00
CoaCAL COAGULATION CALIBRATOR **REF** 210-24-070-00
CoaCPN CONTROL PLASMA NORMAL **REF** 210-24-050-00
CoaCPA CONTROL PLASMA ABNORMAL **REF** 210-24-060-00

Standard laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipette tips.

Plastic test tubes for sample dilution.

Coagulation tubes.

SAMPLES

Collection of samples in accordance with CLSI (NCCLS) (see References 3) using citrated tubes.

ANALYTICAL PROCEDURE CoaDATA SERIES (semi-automated)

Refer to your instrument manual or available application sheets provided for more detailed instrument specific instructions.

Stepwise pipette scheme for manual systems:

50µl Sample, control or standard plasma

50µl CoaPTT

120s Incubation at 37°C (in the measuring block).

50µl CoaCC CALCIUM CHLORIDE as start reagent (pre-warmed at 37°C).

The clotting time in seconds will be recorded automatically by the analyzer. Refer to the instrument operator's manual for further information.

ANALYTICAL PROCEDURE with MANUAL METHOD (see References 4)

- Preheat CoaCC CALCIUM CHLORIDE, liquid at 37°C.
- Pipette into coagulation tube as follows:
 - Plasma 0,050 mL
 - aPTT reagent 0,050 mL
- Incubate for 5 minutes at 37°C.
- CoaCC CALCIUM CHLORIDE, liquid 0,050 mL

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

QUALITY CONTROL

Reliability of test results should be monitored by use of LABiTec Control Plasma CoaCPN and CoaCPA. Control range should be established by the laboratory to determine the allowable variation in day-to-day performance for each level of control. Failure to meet quality control specifications should be investigated and resolved.

CoaCPN CONTROL PLASMA NORMAL **REF** 210-24-050-00

CoaCPA CONTROL PLASMA ABNORMAL **REF** 210-24-060-00

In case the values of Control Plasma are out of the given ranges, check:

Calibration, reagent, stability of sample/controls or analyser.

ATTENTION!

The kit is tested on a Coagulometer CoaLAB 1000 and CoaDATA 4004 (both LABiTec) and ACL 7000 (I.L.). Applications on automatic / semi-automatic or manual coagulometer from other manufacturers may be totally different from what we experienced.

CALCULATION OF RESULTS

For best results, duplicate samples are recommended. APTT results should be reported as clotting time in seconds. Calculate the mean clotting time of duplicate samples and controls. Differences between duplicate results should be less than 5%. Repeat the test if necessary.

REFERENCE VALUES (see References 8)

The results of CoaPTT are influenced by the method of determination of the clot, including reagent and instrument, which may vary from laboratory to laboratory. Generally, the CoaPTT test with normal plasmas originates clotting times between 19 and 37 seconds. However, each laboratory must establish its own normal range by testing individual's representative of the population.

EXPECTED VALUES

Analyser	Method	Normal Values (in seconds)
CoaDATA	Manual	25 - 39
CoaLAB	Automated	25 - 39

ANALYTICAL PERFORMANCES


(validate on CoaLAB 1000, CoaDATA 4004 and ACL 7000 (I.L.))

The performances of the CoaPTT Reagent has been tested with a CoaLAB 1000, CoaDATA 4004 and ACL 7000 Coagulometer. The data, while representing the characteristics of the product, could be different for each laboratory and for different coagulometers.

Heparin Therapy: (see References 5,6)

Therapeutic range: 0,2-0,4U/ml.

We prepared solutions of normal plasma with scalar amounts of unfractionated Heparin to verify the sensitivity of the Reagent. The experimental results obtained in our laboratory with ACL 7000 (I.L.), with specific lot of aPTT, are shown in this graph.

The results reported here are indicative and representative of the characteristics of the reagent on our instruments, but in no case should be used as such in the individual laboratories. Since there are many variables that can affect results (eg. use of different heparins), each laboratory should establish its own therapeutic reference range.

Intrinsic Factor Sensitivity: The experimental tests have shown that the CoaPTT Reagent is sensitive to changes in the factors of the intrinsic pathway.

Method Limitations: (see Reference 2)

Interferences: (see References 2)

Interference test criterion: recovery ± 15% of initial value.

No interference found on samples with:

- total bilirubin up to 40 mg/dL;
- hemoglobin up to 600 mg/dL;
- lipemia [Intralipid ®] up to 1000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples. The results obtained are following:

	CoaLAB	CoaDATA			ACL 7000	
Sample	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%
Human Plasma 1	28,2 ± 0,55	1,04	38,57 ± 1,21	1,0	29,5 ± 0,60	1,1
Human Plasma 2	37,72 ± 0,6	1,09	60,75 ± 2,38	1,1	37,9 ± 0,96	1,2

Run-to-run Precision: determined for 5 days with 20 replications for each day, for 2 samples.

The results obtained are the following:

	CoaLAB	CoaDATA			ACL 7000	
Sample	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%	Average (seconds) ± 2s	CV%
Human Plasma 1	42,1 ± 0,93	1,11	43,6 ± 1,05	1,2	38,4 ± 0,87	1,11
Human Plasma 2	57,3 ± 1,13	0,99	58,2 ± 0,82	0,703	49,5 ± 1,19	1,2

Accuracy: a group of 20 plasma has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression: $y = 1,0133x - 0,900$

Correlation coefficient: $r = 0,9953$ $n = 20$

WARRANTY

This product is warranted to perform in accordance with its labeling and literature. LABiTec disclaims any implied warranty of merchantability or fitness for any other purpose. Purchaser must calibrate and determine the suitability of LABiTec's products for their specific applications. In no event will LABiTec be liable for any consequential damages arising out of aforesaid express warranty.

ANALYTICAL PROCEDURE CoaLAB SERIES (automated Instrument)

Refer to CoaLAB instrument user manual

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